

# *Increase your LC resolution*

*Large capacity column oven  
(up to 20 columns)*

*Reaching 200,000 up to  
500,000 plates*

*Conventional columns  
and LC equipment*

*Dedicated preheater*

*Dedicated postcooler*

***Polaratherm™  
Series 9010  
HIGH RESOLUTION***



The Polaratherm Series 9010 High Resolution enables high efficiency LC separations. This is achieved on standard LC equipment by coupling conventional LC columns and performing the separation at elevated temperature.

**SANDRA SELERITY TECHNOLOGIES INC.**

USA, Canada, and Latin America: [www.selerity.com](http://www.selerity.com) • [jodyclark@selerity.com](mailto:jodyclark@selerity.com) - Other countries: [www.richrom.com](http://www.richrom.com) • [sales@richrom.com](mailto:sales@richrom.com)